

## DNA isolation from mature endosperm

*(Based on DNA isolation protocol by Mauricio Lopes)*

*The conditions for DNA isolation from mature endosperm should be optimized to avoid precipitation of starch with genomic DNA. To overcome this problem, a sarcosyl buffer is used and samples are not heated during extraction. To improve intact DNA yields, kernels (30 DAP or older) can be stored at – 80 C or lyophilized, which also facilitates the grinding and flour weighting.*

### Procedure

Dissect the endosperm removing the pericarp. Embryonic tissues can also be used;  
Lyophilize the tissue – OPTIONAL- for 24 hours (6 hours @ -20 C and 18 hours @ 0 C);  
Move it to a pre-chilled mortar and grind it making a fine flour;  
Weigh 50-200 mg of flour in a eppendorf tube;  
Add 600 ul of extraction buffer (see below);  
Vortex, mixing up the flour and buffer well;  
Let it stand @ RT for 10 min.;  
Add 600 ul of Phenol-Chlorophorm, vortex it vigorously;  
Transfer to a PHASE LOCK tube, if available;  
Spin @ 14000 rpm for 10 min.;  
Transfer supernatant to a new tube, add 5 ul of RNase;  
Incubate @ 37 C for one hour;  
Add 450 ul of cold isopropanol;  
Mix and incubate @ -80 C for 30 min.;  
Spin @ 14000 rpm for 10 min.;  
Pour off supernatant and wash the pellet with 500 ul of 70% Ethanol;  
Dry the pellet @ RT;  
Add 50 ul of water, resuspend the pellet overnight @ 4 C;  
Spin for 15 sec.  
Run a gel with Lambda DNA standards to check quality and determine concentration.

### Buffer and reagents

<b>Extraction Buffer</b>	100 mM TRIS pH 8.5 100 mM NaCl 200 mM EDTA pH 8.0 1% Sarcosyl
<b>RNase</b>	10 mg/ml (Boil to inactivate any DNase activity)
<b>Isopropanol</b>	Kept @ -20 C